

KERN ORT1RS

Area of application: Industry/Pharmaceuticals/Laboratory



Measuring system

Accuracy:	0.0002 nD 0.1 %
Division:	0.0001 nD 0.1 %
Measuring range:	1.3000 nD - 1.7000 nD 0.1 % - 95 %
Scales:	Refractive index Brix

Packaging & shipping

Delivery:	24 h
Dimensions packaging (WxDxH):	380 x 255 x 165 mm
Gross weight:	3 kg
Net weight:	1,950 kg
Shipping method:	Parcel service